

Assay Data Review and Release Work Instruction

Loblolly Ridge Therapeutics | Discovery Biology | 2025-10-20

Work instruction WI-DB-007, revision B. It describes how an in vitro assay run is completed, checked, uploaded, reviewed and released, from the instrument to the assay registry and the R&D Data Hub. It applies to all Discovery Biology and DMPK assay runs at the Woodlands laboratory.

1. Purpose and principles

Assay results drive compound design and project decisions, so every released result must be traceable to its run, its raw data and the people who produced and reviewed it. Research at the Woodlands laboratory is non-GxP, but we apply attributable, contemporaneous and reviewed record-keeping as good practice. FDA describes data integrity as the completeness, consistency and accuracy of data; audit trails in the ELN and the registry let changes be reconstructed.

2. Where the data go

Results flow from the instrument and chromatography data system into the ELN experiment, then to the assay registry, and nightly into the R&D Data Hub, where they are joined to compound identifiers from the compound database. Only released results reach the Data Hub dashboards.

3. Roles

Scientist: runs the assay, records the ELN entry, applies acceptance criteria and uploads the run.

Reviewer: a second qualified scientist who did not perform the run; signs off or returns it.

Research Informatics: maintains the registry and reviews correction requests.

4. Steps

4.1 Complete the run and the ELN entry on the same or the next working day. Attach the instrument export and plate map, and record reagent lots and deviations from the method.

4.2 Apply the acceptance criteria for the assay (controls, reference compound, plate statistics). Criteria for each assay are in its method record.

4.3 If the run fails, record why. Repeat the plates when material allows and record the outcome as fail-repeat, or record fail-invalid if the run cannot be repeated.

4.4 Export the run from the ELN and import it into the assay registry. This is a manual step done by the scientist; check that the plate count and compound IDs match the ELN entry.

4.5 Request review in the registry. The reviewer checks the raw data, the calculation and the acceptance criteria, then signs off or returns the run with comments.

4.6 On sign-off, release the run. If a run is invalidated, choose a reason code from the registry pick list and add a comment.

5. Corrections after release

A released result cannot be edited. To correct it, submit a registry correction request with the reason and the evidence; Research Informatics reviews it and the reviewer of the original run confirms the change. The registry keeps the original and corrected values.

6. Timeliness

Aim to upload a run within two working days of the ELN entry and to complete review within three working days of the request. Project teams read the Data Hub dashboards each morning, so an unreleased run is invisible to them.

7. Revision history

Rev. B, October 2025 – fail-repeat and fail-invalid outcomes defined; reviewer independence stated; correction requests routed to Research Informatics.

Rev. A, March 2025 – first issue with R&D Data Hub phase 1.